

Organo Corporation

The Power of Purification

Driving Industrial Innovation with Advanced Water Systems

US\$1bIn

Annual Sales

75+

Years in Business

2,500+

Employees

8

Countries

We are Organo, Japan's leading provider of water treatment and purification systems. Our engineering expertise is at the forefront of purification through systems that deliver water of impeccable purity, and much more. Backed by leading-edge R&D, we deliver UPW and other solutions that elevate quality, create a competitive edge, and meet the evolving customer needs of various industries.

Pharma Business

Serving the Pharmaceutical Industry for 40 Years

Top market share in Japan

Over 550+ systems delivered

Organo has established itself as a leader in pharmaceutical water treatment, developing innovative solutions using advanced technologies. The company provides comprehensive pharmaceutical water equipment and consulting services to meet evolving industry regulations for pharmaceutical, cosmetic, and medical device sectors in Japan and abroad.

Designed and Engineered in Japan, Assembled and Delivered Globally

Organo supplier networks cultivated in Asia and Europe form our competitive advantage. Our technology, quality, and services satisfy customers worldwide.



System Engineering and Design

Tokyo-based strategic command center

Organo headquarters in Tokyo is the main base for our comprehensive range of engineering development and project planning services.


Vietnam Factory

Solutions from APAC to the world

As a pivotal link in our global network, the Vietnam factory has been delivering tailor-made solutions to the Asia-Pacific region and beyond since 2010. It stands as a crucial hub for the expansion of our global presence in the pharmaceutical sector.


Product Information

PWG Series - Water Purification System



Model	PWG 1-0.5	PWG 1-1	PWG 1-2	PWG 1-3
Production Capacity	0.55m ³ /h	1.1m ³ /h	2.2m ³ /h	3.3m ³ /h
Dimensions (mm)	W2,100 x D1,950 x H1,970	W2,100 x D1,950 x H1,970	W2,100 x D1,950 x H1,970	W2,720 x D2,600 x H2,100
Operating Weight	2.5 tons	2.5 tons	3.4 tons	4.0 tons

PWD Series - Purified Water Storage and Distribution System



Model	PWD1-030	PWD1-050	PWD1-070	PWD1-110
Supply Capacity	~5m ³ /h	10-15m ³ /h	15-30m ³ /h	30m ³ /h-
Dimensions (mm)	W1,600 x D800 x H2,100	W1,600 x D800 x H2,100	W1,600 x D800 x H2,100	W1,600 x D800 x H2,100
Operating Weight	0.8 ton	0.8 ton	0.9 ton	1.0 ton

Contact Us

Yuichiro Toyoda | Senju Hara

pharma-usa@organo.co.jp

Visit us at www.organo.co.jp/english/business/medical/